


Baxter

GLOBAL TRAVEL CLUB GUIDE



Your Destination
for Traveling with
Peritoneal Dialysis





Introduction	4
The Benefits of Global Travel Club	6
Useful Contacts	7
Travel Order Process	8
Fee Schedule	10
Summary of Responsibilities	12
Cruising: PD at Sea.....	14
Travel Considerations:	
» Continuous Ambulatory Peritoneal Dialysis (CAPD)	16
» Automated Peritoneal Dialysis (APD)	17
Traveling With Your Cyclor:	
» Homechoice Cyclor	20
» Amia/Kaguya Cyclor	21
Documentation for Travel Orders	22
Other Considerations To Keep in Mind	24
The PD Travel Checklist	25
Frequently Asked Questions.....	26
Useful Contacts	29
Appendix:	
» Connector Guide	30
» Customize Additional Information	31

Whether it's to visit family, spend time in a second home or cottage, explore a foreign land, or just a case of wanderlust, the freedom and flexibility of Peritoneal Dialysis (PD) opens up a world of possibilities. Over the years, Baxter has helped thousands of people fulfill their travel dreams, whether it be a trip for the weekend or the trip of a lifetime.

At Baxter, we understand that traveling while on dialysis comes with a unique set of challenges. That's why we support our patients with the Global Travel Club.

- » To help arrange delivery of your PD supplies directly to wherever your travels take you
- » To help guide you through the process, from start to finish
- » And perhaps most importantly, to help you live a life well-traveled



Our experienced team is ready to assist you, so please give us a call if you have any questions about this service.

READY  YOUR JOURNEY AWAITS!





THE BENEFITS OF GLOBAL TRAVEL CLUB

One of the many great benefits of choosing PD is that it gives you freedom to travel almost anywhere in the world. Through Baxter's global contacts and partnerships, we're able to safely ship your PD supplies directly to your destination, simplifying and streamlining the process.

When factoring in how long your trip is and the amount of PD supplies required, it's easy to see just how valuable the Global Travel Club can be.

LITERS PER DAY

LENGTH OF TRIP

TOTAL LITERS NEEDED

12*



3
DAYS



36



12*



7
DAYS



84



12*



14
DAYS



168



*12 liters per day is a sample prescription used for illustrative purposes only. Your specific numbers may vary.



USEFUL CONTACTS

You're bound to have questions. Not to worry, we've got answers. Here's a handy list of resources you may find helpful.

GLOBAL TRAVEL CLUB: 1-800-XXX-XXXX

Phone Number 1: 1-800-XXX-XXXX

Phone Number 2: 1-800-XXX-XXXX

Phone Number 3: 1-800-XXX-XXXX

Phone Number 4: 1-800-XXX-XXXX

TRAVEL ORDER PROCESS

There's nothing like the excitement of planning an upcoming trip. And part of that preparation includes securing your PD supplies so you can travel with confidence. To help ensure you don't experience any unexpected turbulence along the way, our Global Travel Club is here to make your travels go smoothly.

Q

Ok, that sounds good... but what exactly do I have to do? And how? And by when?

A

All great questions! We'll break it down in more detail, step-by-step. And as always, if you need help, please don't hesitate to contact a Global Travel Club Specialist at any time.



STEP 1: BEFORE MAKING DEFINITE TRAVEL ARRANGEMENTS AT LEAST 75 DAYS IN ADVANCE

- » Talk to your clinical PD team (your PD nurse or doctor) about your potential travel plans
- » Get the necessary permissions from your clinical PD team to travel to your destination for a specified length of time
- » **Your travel order must be approved by your clinical PD team.** Baxter will confirm with you that the appropriate approvals are obtained before collecting any shipping and handling fees, if applicable



STEP 2: CONTACT THE BAXTER GLOBAL TRAVEL CLUB A MINIMUM 60 DAYS' ADVANCE NOTICE IS REQUIRED (THIS MAY VARY BY COUNTRY)

- » Tell us about your destination so your Global Travel Club Specialist can confirm that Baxter is able to deliver your required supplies
- » Provide travel dates, destination address and contact information as well as any special instructions
- » **Before** booking your travel, you will need to secure signatures on two forms:
 - Travel waiver — to be signed by both patient and kidney clinic
 - Prescription form — to be approved by your physician
- » Additional information may be required (passports, affidavits, medical history report, customs paperwork, etc.)
- » Your Global Travel Club Specialist will provide a quote for any costs, if applicable



STEP 3: CONTACT YOUR DESTINATION


Make sure your destination, whether it be a hotel, cruise ship or private residence, can receive and properly store your supplies prior to your arrival

- » Confirm that the temperature and cleanliness of the storage room at your destination are appropriate
- » Baxter will try to deliver supplies immediately prior to your arrival date



STEP 4: PLACE THE ORDER

- » Baxter will confirm your order has been received
- » As the patient, you are responsible for verifying that the details of the order are correct
- » **Once you receive that confirmation and it's correct, then you can book your travel plans**
- » Baxter will collect payment for any costs, should they apply



STEP 5: PREPARE TO TRAVEL WITH YOUR CYCLER

- » If you're planning on traveling with your cycler, make sure you have the required power adapter and/or transformer depending on your specific cycler
- » Consult your clinic as to the best way to transport your cycler — they may have a special carrying case you should use
- » **Homechoice Cycler or Homechoice Pro Cycler:** requires a transformer plus a power adapter
- » **Homechoice Claria Cycler or Amia/Kaguya Cycler:** does not require a transformer and only needs a power adapter



STEP 6: RIGHT BEFORE YOU TRAVEL

- » Baxter will provide an expected delivery date of when your product will arrive at the destination
- » Before departing for your trip, please check with your destination that delivery has been made



STEP 7: LAST DAY OF YOUR TRIP

- » When leaving your destination, it is your responsibility to dispose of waste or unused product; Baxter is unable to pick up any unused products
- » Please note: Baxter is unable to provide any refunds for unused product



FEE SCHEDULE

Shipping medical supplies, especially internationally, is increasingly complex in today's constantly changing environment. Depending on the services required and your destination, some fees, such as freight, customs, and handling charges, may be associated with the delivery of your product. These fees will be your personal responsibility. Your Global Travel Club Specialist will be able to provide a quote specific to your travel itinerary, if fees are applicable to you.

What Will It Cost?

DOMESTIC

For travel within the country, Baxter will cover the product freight costs associated with your travel

INTERNATIONAL

(Baxter Network)

For international destinations, freight charges are kept to a minimum because of our network of Baxter locations worldwide. Your goods are shipped from the closest possible Baxter or distributor location to the travel destination

INTERNATIONAL

(Outside Baxter Network)

In some instances, Baxter will be unable to source the product locally. In these cases, your PD supplies will need to be shipped directly from your home country at an additional cost, which is your responsibility

CRUISES

Due to the complexities associated with the delivery of PD supplies to cruises, an additional fee will apply

EXPEDITED

If you are unable to provide the required minimum 60-day lead time, an expedited fee will apply





→ **PRICE** ←

→ **PRICE** ←

→ **PRICE** ←

→ **PRICE** ←

→ **PRICE** ←

Please Note: Additional fees, such as customs clearance fees, payment of any local duties or taxes on PD fluids, consumables, storage fees, returns and proper disposal of extra product, may arise on a case-by-case basis. As these fees are not Baxter-associated, **they are solely the responsibility of the patient traveling.** Baxter is unable to advise on what these costs may be; patients must be prepared to pay these fees as required on site.



SUMMARY

OF RESPONSIBILITIES

To place an order, you'll work with your clinical team as well as our experienced Global Travel Club Specialists, who will help you navigate the process.

(YOUR) PATIENT'S RESPONSIBILITIES

BEFORE making definite travel plans, talk to your clinical PD team to obtain the necessary permissions to travel

Contact the Global Travel Club, **AT LEAST 60 DAYS IN ADVANCE**, to begin the ordering process (lead time varies by country)

Confirm with your travel destination that they can properly store supplies prior to your arrival

Sign travel waiver and return to Baxter

Pay any applicable fees

CLINICIAN'S RESPONSIBILITIES

Provides necessary permissions for you to travel

Trains you on PD procedures, if necessary

Provides a carrying case, if available, for you to transport your APD cyclers

GLOBAL TRAVEL CLUB'S RESPONSIBILITIES

Coordinates delivery of PD supplies to your destination

Provides you with a quote and then collects payment of Baxter-related delivery fees, if applicable

Provides you with an expected delivery date of when your PD supplies will arrive at your destination

IMPORTANT! Please contact Baxter to begin making travel arrangements as early as possible. The more time we have to plan, the smoother the travel arrangements can go, and this is also the best opportunity to minimize expenses. The sooner you contact Baxter, the better, and make sure it's **BEFORE** you confirm your travel plans.



C R U I S I N G

PD AT SEA



The Global Travel Club is here to help make cruising with PD smooth sailing. While we strive to make this process as easy as possible, there are logistical challenges surrounding the coordination of your PD supplies with respect to cruising. Therefore, an additional fee will apply. Please make sure to have the following information readily available for the Global Travel Club Specialist at the time the order is placed:

- 1 BOOKING NUMBER
- 2 CRUISE LINE
- 3 SHIP NAME
- 4 CABIN/STATEROOM NUMBER
- 5 SAIL DATE AND TIME
- 6 RETURN DATE AND TIME
- 7 ADDRESS OF PORT OF DEPARTURE/ADDRESS OF WAREHOUSE REQUIRED FOR DELIVERY
- 8 PORT AGENT/SHIPPING SERVICES CONTACT INFORMATION

IMPORTANT: A MINIMUM 60 DAYS' NOTICE IS REQUIRED FOR ALL CRUISE ORDERS.

Please check that the port agent is able to accept deliveries on the required day(s). If the product is then kept in storage and charged, you will be responsible for those additional charges. These charges are assessed by the cruise line or port authority.

TRAVEL

CONSIDERATIONS

CAPD

When traveling as a CAPD patient, which means you perform manual exchanges during the day, please take note of the following important information.

- » In order to travel while on CAPD therapy, you'll need solution bags and other disposable products such as **MiniCap** disconnect caps, clamps, etc.
- » Remember to either arrange for or personally bring all of the disposable products needed
- » In many instances, it may be easier for you to bring disposable products vs. having them delivered, especially for shorter trips



APD

If you're an APD patient, which means you do PD therapy with a cycler at night, there's some important information you'll want to be aware of when traveling.

- » In order to travel while on APD therapy, you'll need a transformer (in some cases), a power adapter, solution bags, the disposable set, and other disposable products such as **MiniCap** disconnect caps, etc.
- » Remember to either arrange for or personally bring all of the disposable products needed
- » In some instances, it may be easier for you to bring disposable products vs. having them delivered, especially for shorter trips
- » You'll need to be aware of the solution bag connection differences (Spike or Luer) between countries (see page 30 for more information)
- » You may want to consider whether or not back-up CAPD supplies should also be delivered to your destination, in the event of APD cycler issues — in many situations, this is highly recommended



TRAVELING WITH YOUR CYCLER



If you're planning to travel with your cycler, you might be wondering, "Wait, what do I need to bring with me for my specific cycler?"

There are a few different cycler options: the **Homechoice** Cycler, the **Homechoice Pro** Cycler, the **Homechoice Claria** Cycler, and the **Amia/Kaguya** Cycler. Depending on which one you use and if you're traveling internationally, a transformer and/or power adapter may be required, which you must supply.

- » **Homechoice** cycler or **Homechoice Pro** cycler: requires a transformer plus a power adapter
- » **Homechoice Claria** cycler or **Amia/Kaguya** Cycler: does not require a transformer and only needs a power adapter

Please note: It is your responsibility to transport your cycler to your destination as Baxter does not ship these machines. Please consult with your clinical team to determine the best way to transport your cycler.





Traveling With Your Cyclor

HOMECHOICE

There are some considerations to be mindful of when planning to travel with a **Homechoice** Cyclor.

- » When traveling internationally with a **Homechoice** cyclor or **Homechoice Pro** cyclor, **a transformer plus power adapter is required.**

Your **Homechoice** cyclor voltage: **110V** **220V** (Check one)

Your **destination's** required voltage: **110V** **220V** (Check one)

If these two numbers are not the same, and you use a **Homechoice** or **Homechoice Pro** cyclor, you will require a transformer.

- » **Homechoice Claria** cyclor does not require a transformer and only needs a power adapter for international travel
- » You should check with your airline for possible carry-on restrictions, as these vary by airline
- » Please see the "Device Travel Letter" in the "Documentation for Travel Orders" section for a sample of the letter you can provide to the airlines



HOMECHOICE OR HOMECHOICE PRO CYCLER

HOMECHOICE CLARIA CYCLER

US IMPERIAL

METRIC

US IMPERIAL

METRIC

HEIGHT

7 in

17.78 cm

7.6 in

19.4 cm

WIDTH

19.5 in

49.53 cm

18.4 in

46.7 cm

DEPTH

15.7 in

39.9 cm

15.2 in

38.7 cm

WEIGHT

27 lb

12.3 kg

29.8 lb

13.5 kg

AMIA/KAGUYA

If you are using an **Amia/Kaguya** Cycler, there are additional travel considerations to be aware of.

- » The **Amia** cycler is only available in Canada and the U.S., and the **Kaguya** cycler for PD, only in Japan. Therefore, if you are traveling with your **Amia/Kaguya** cycler outside these countries, there is no technical support available
- » Some countries carry spike solution bags that are not compatible with an **Amia/Kaguya** cycler; if you plan to travel to one of these areas, contact your clinical team ahead of time
- » When traveling internationally with an **Amia/Kaguya** cycler, a transformer is not required; only a power adapter is needed
- » You should check with your airline for possible carry-on restrictions, as these vary by airline
- » Please see the “Device Travel Letter” in the “Documentation for Travel Orders” section for a sample of the letter you can provide to the airlines



AMIA/KAGUYA CYCLER

US IMPERIAL	METRIC
7.20 in	18.3 cm
14.96 in	38 cm
11.42 in	29 cm
18.08 lb	8.2 kg



DOCUMENTATION FOR TRAVEL ORDERS

TRAVEL ORDER FORM

- » Your clinical PD team will work with you to fill out and submit the required Travel Order Form

TRAVEL WAIVER

- » Baxter has a travel waiver that you will be required to complete, sign, and return to Baxter. Your Global Travel Club Specialist will share this with you

DEVICE TRAVEL LETTER

- » To help avoid airport delays, you can present this letter to the airlines, which contains information about your cyclor. This letter and other helpful travel documents are available on our website at GlobalTravelClub.Baxter.com

Travel Order Form



Device Travel Letter

DOCTOR'S LETTER

- » This is an important letter as it outlines some of the specifics of your PD therapy, as well as your trip details. Use this letter (below) to help place your Global Travel Club order, to present to Customs and Immigration (if traveling internationally), and to present to the airlines to help avoid airport delays



(TO BE PRINTED ON HOSPITAL LETTERHEAD)

PATIENT NAME

PATIENT ADDRESS

RE: **(INSERT PATIENT NAME)** of **(INSERT PATIENT ADDRESS)**

To whom it may concern,

(INSERT PATIENT NAME) is traveling to **(COUNTRY OF TRAVEL)** on **(INSERT DATES OF TRAVEL)**.

(INSERT PATIENT NAME) is a kidney failure patient under the care of this medical facility and uses Peritoneal Dialysis as a life support system. The following supplies are required for his/her treatment.

- » **(INSERT PRODUCT LIST AS ADVISED BY BAXTER TRAVEL SPECIALIST)**

These products have no commercial value and are for the sole use of the patient named in this letter.

Signed,

- » **(INSERT DOCTOR'S SIGNATURE AND DOCTOR'S CONTACT INFORMATION)**



OTHER CONSIDERATIONS TO KEEP IN MIND



PRODUCT

- » Products may look slightly different from country to country; talk to your PD clinical team to discuss if you need additional training before your trip. The product you receive at your travel destination may have different product code numbers on the box, and depending on your destination, may be in a different language
 - It's important that you familiarize yourself with the destination packaging, to make sure you're using the correct product(s)
- » In some cases, products that Baxter offers are not available in all countries. It's possible that you may be required to bring your own disposable products; your Global Travel Club Specialist will tell you if this applies to you. You may also have the option of shipping the product from your home country for an additional charge
- » Baxter will not pick up unused product from a destination, nor will we provide refunds


➔ You are responsible for proper disposal of unused product ➔



DELIVERY

- » Baxter will deliver a maximum supply of four weeks at a time
 - If you are staying longer, you may need to arrange multiple deliveries
 - Your Global Travel Club Specialist will give you details specific to your situation
- » **In some situations**, Baxter can only deliver to the local airport, and you will be responsible for transporting your product from the airport to your local accommodations
- » It is recommended that you contact your country's consulate to confirm customs procedures
- » Remember to bring dressings and medication supplies with you
- » Familiarize yourself with voltage and transformer requirements
- » Be prepared for "unexpected" expenses
- » If you are an APD patient, it is strongly suggested that you also bring some CAPD supplies as back-up
- » In some instances, we will be unable to meet your international travel needs. **Please contact Baxter early in your travel planning and before you finalize your plans to avoid any unexpected surprises**

THE PD TRAVEL CHECKLIST



In addition to the regular items that you usually take (e.g., clothes, toiletries, books) here is a checklist to help you with your PD packing:

TO DO BEFORE LEAVING:

- Confirm with your destination that your supplies have arrived
- Gather a list of contact numbers: Baxter, your clinician, etc.
- Understand if your supplies will be delivered directly to your lodging destination, or if they need to be picked up by you at the airport

SUPPLIES TO TAKE WITH YOU:

- PD cyclor, if needed
- Transformers and power adapters for international travel, if needed
- Medications: Non-Baxter medications needed for your trip should be in their original packaging
- PD supplies, dressings and medication supplies: Please ensure you take all small items you will require for your treatment, as these are sometimes not shipped; review your supply list to be sure
- Amia** Cassettes and Supplies: You may need to take your **Amia** cassettes and supplies with you depending on the destination. Check with your Global Travel Club Specialist

DOCUMENTS TO TAKE WITH YOU:

- You **MUST** carry your “Doctor’s Letter” outlining your condition, in order to travel. This letter will also be required at Customs and Immigration, if you plan on international travel
- Copy of prescriptions
- Baxter’s order confirmation
- Homechoice** Cyclor or **Amia** Cyclor “Device Travel Letter”



FAQ

FREQUENTLY ASKED QUESTIONS

»» WHAT KIND OF LEAD TIME IS REQUIRED?

In order to properly arrange for the delivery of your PD supplies when traveling, the Global Travel Club requires a minimum of 60 days' advance notice.

»» WHY DO YOU NEED ADVANCE NOTICE?

The more notice we have, the better we can assist you in your travel plans. Shipping medical supplies, especially internationally, is a complex process. By knowing your plans in advance, we can advise you:

- »» If the destination is supported by the Baxter Global Travel Club
- »» What, if any, additional costs may be required so that you can factor these into your plans
- »» What paperwork is needed for travel approval by your clinical team
- »» What paperwork is required by your destination

»» WHY DO I NEED TO SIGN A TRAVEL WAIVER?

Planning travel has some risks involved, including canceled or changed plans. The travel waiver exists to make sure that you understand the risks and are willing to assume the liability involved with planning travel.

»» WHAT IF MY TRAVEL PLANS CHANGE?

Sometimes plans change. We understand. Please contact a Global Travel Club Specialist as soon as possible so we can assist you.



»» HOW MUCH WILL THIS COST?

We know that cost is an important factor when planning a trip. Depending on the services required and your destination, some fees, such as freight, customs and handling charges may be associated with the delivery of your product. Your Global Travel Club Specialist will be able to advise you on any costs, if they are applicable to you.

»» WHAT IF I ARRIVE AND MY SUPPLIES ARE NOT THERE?

While we take every measure to ensure your supplies arrive on time, sometimes things do happen that are out of our control. If you are advised that your supplies have not arrived or you can't obtain confirmation of delivery, contact a Global Travel Club Specialist immediately.

»» HOW DO I DISPOSE OF UNUSED PRODUCT?

Unused disposable products like APD sets, caps, etc. should be disposed of in accordance with local laws. Unused solution bags should be drained first, then the containers should be disposed of in accordance with local laws.

»» WHOM DO I CONTACT IF I HAVE QUESTIONS?

Having questions is all a part of the process. And we're happy to help. That's why we're here. Our dedicated Global Travel Club Specialists are ready to answer whatever's on your mind.

GLOBAL TRAVEL CLUB: 1-800-XXX-XXXX







USEFUL CONTACTS

As you made your way through this guide, we hope it provided the information you need to travel with confidence. However, we understand there are bound to be questions. The list below contains the contact information of some resources you may find useful. And as always, our experienced Global Travel Club Team is here to help you navigate this process.

GLOBAL TRAVEL CLUB: 1-800-XXX-XXXX

Phone Number 1: 1-800-XXX-XXXX

Phone Number 2: 1-800-XXX-XXXX

Phone Number 3: 1-800-XXX-XXXX

Phone Number 4: 1-800-XXX-XXXX

YOUR PD NURSE'S CONTACT INFORMATION



SAFE TRAVELS!

Every trip takes you to new places. Whether it's exploring a new city or just taking the time to relax, traveling feeds your mind, body and soul. And Baxter's Global Travel Club can help you get there.

SO, GO AHEAD — THE WORLD IS CALLING, AND YOUR ADVENTURE AWAITS!

For more information about traveling while on PD, please visit GlobalTravelClub.Baxter.com.

Connector types vary by country.
Please talk to your clinician to determine your needs.



Countries That Offer Spike Connectology

Asia Pacific	Europe, Middle East and Africa	Americas		
CHINA	JORDAN	ARGENTINA	ECUADOR	NICARAGUA
JAPAN	SOUTH AFRICA	BELIZE	EL SALVADOR	PANAMA
		BRAZIL	GUATEMALA	PARAGUAY
		CHILE	HONDURAS	URUGUAY
		COLOMBIA	MEXICO	VENEZUELA

Countries That Offer Luer Connectology

Asia Pacific	Europe, Middle East and Africa		Americas
AUSTRALIA	ALBANIA	LITHUANIA	BARBADOS
BANGLADESH	ALGERIA	LUXEMBOURG	BERMUDA
BRUNEI	AZERBAIJAN	MACEDONIA	CANADA
CHINA	AUSTRIA	MALTA	COSTA RICA
HONG KONG	BAHRAIN	MONTENEGRO	DOMINICAN REPUBLIC
INDIA	BELGIUM	MOROCCO	ECUADOR
INDONESIA	BOSNIA AND HERZEGOVINA	NETHERLANDS	GUYANA
JAPAN	BULGARIA	NORWAY	JAMAICA
MALAYSIA	CROATIA	OMAN	PERU
MYANMAR	CZECH REPUBLIC	POLAND	PUERTO RICO
NEPAL	DENMARK	PORTUGAL	TRINIDAD AND TOBAGO
NEW CALEDONIA	ESTONIA	QATAR	UNITED STATES
NEW ZEALAND	FINLAND	ROMANIA	VENEZUELA
PHILIPPINES	FRANCE	RUSSIA	
SINGAPORE	GEORGIA	SAUDI ARABIA	
SOUTH KOREA	GERMANY	SERBIA	
SRI LANKA	GREECE	SLOVAKIA	
TAIWAN	HUNGARY	SLOVENIA	
THAILAND	IRAN	SOUTH AFRICA	
VIETNAM	IRAQ	SPAIN	
	IRELAND	SUDAN	
	ISRAEL	SWEDEN	
	ITALY	SWITZERLAND	
	JORDAN	TUNISIA	
	KAZAKHSTAN	TURKEY	
	KUWAIT	UKRAINE	
	LATVIA	UNITED ARAB EMIRATES	
	LEBANON	UNITED KINGDOM	
	LIBYA		



COUNTRY TO CUSTOMIZE WITH PI, ADDITIONAL INFORMATION, ETC.

Traveling With Your Cyclor

HOMECHOICE

There are some considerations to be mindful of when planning to travel with a **Homechoice** Cyclor.

- » When traveling internationally with a **Homechoice** cyclor or **Homechoice Pro** cyclor, **a transformer plus power adapter is required.**

Your **Homechoice** cyclor voltage: **110V** **220V** [Check one]

Your **destination's** required voltage: **110V** **220V** [Check one]

If these two numbers are not the same, and you use a **Homechoice** or **Homechoice Pro** cyclor, you will require a transformer.

- » **Homechoice Claria** cyclor does not require a transformer and only needs a power adapter for international travel
- » You should check with your airline for possible carry-on restrictions, as these vary by airline
- » Please see the "Device Travel Letter" in the "Documentation for Travel Orders" section for a sample of the letter you can provide to the airlines



HOMECHOICE OR HOMECHOICE PRO CYCLOR

HOMECHOICE CLARIA CYCLOR

US IMPERIAL

METRIC

US IMPERIAL

METRIC

HEIGHT

7 in

17.78 cm

7.6 in

19.4 cm

WIDTH

19.5 in

49.53 cm

18.4 in

46.7 cm

DEPTH

15.7 in

39.9 cm

15.2 in

38.7 cm

WEIGHT

27 lb

12.3 kg

29.8 lb

13.5 kg





Baxter

Baxter, Amia, Homechoice, Homechoice Claria, Homechoice Pro, Kaguya and MiniCap are trademarks of Baxter International Inc. Any other trademarks, product brands or images appearing herein are the property of their respective owners.

GLBL/MG2/20-0020 09/20